
Benchtop qNMR Spectroscopy **A Journey Across** **Three Continents**

SNUG PG Meeting
January 12th, 2023

Alexander Maier, PhD



Company Snapshot



Sell Direct in USA, Canada, Germany, and other parts of Europe;
Distributor Network of 45 dealers in the UK and the rest of world.



Technology Leader

- Products in LF, HF and MRI
- Patent protected technology



MR Research Leader

- Exciting international research collaborations with government & industry



Growing Global Business

- > 200 employees
- 1000+ customers worldwide

Nanalysis Scientific Corp.

TSXV : NSCI | OTC : NSCIF | FRA : 1N1

A global product development and manufacturing company designing a suite of leading-edge, high-performance magnetic resonance technologies to provide innovative solutions to customers in academia, government and industry.



Miniaturization: *accessible, affordable, automatable*

Magnet platform

- NdFeB magnets arranged in patented, highly homogenous compact Halbach based arrays
- Cryogen free

Electronics

- Condensed to a small circuit board stack
- Embedded computer for data acquisition and processing



Nanalysis Benchtop NMR Spectrometers



60e



60PRO



100e



100PRO

	60 MHz (1.4 T)		100 MHz (2.35 T)	
Field strength	60 MHz (1.4 T)		100 MHz (2.35 T)	
Magnet	NdFeB		NdFeB	
Weight	26 kg		110 kg	
Cryogen Free	Y		Y	
Shielding	<2 Gauss outside enclosure		<2 Gauss outside enclosure	
Operating Temperature	18-26°C		18-26°C	
Touchscreen interface	Y		Y	
Nuclei	$^1\text{H}/^{19}\text{F}$	$^1\text{H}/^{19}\text{F}/\text{X}$ X: ^{11}B , ^{13}C , ^{31}P , ...	^1H	$^1\text{H}/\text{X}$ X: ^{11}B , ^{13}C , ^{19}F , ^{31}P , ...
Experiments	1D, T_1 , T_2 , COSY, JRES	1D, T_1 , T_2 , COSY, JRES DEPT, HSQC, HETCOR, HMBC, Decoupling	1D, T_1 , T_2 , COSY, JRES	1D, T_1 , T_2 , COSY, JRES DEPT, HSQC, HETCOR, HMBC, Decoupling

Geographic Distribution of NMR



	Country	%
1	United States	32
2	Germany	11
3	China	10
4	Japan	10
5	UK	6
6	France	4
7	India	4
8	Switzerland	3
	Total	80

Quantitative NMR (qNMR)

NMR is inherently quantitative

Integration area of each signal is directly proportional to the number of chemically equivalent nuclei giving rise to that signal

Components do not need to be separated

Provides a snapshot of everything in solution

A single product resonance can be used to quantify


The analysis is non-destructive

The analyte can be recovered

Calibrant methods:

- 1) Internal calibrant
- 2) External calibrant
- 3) Calibration curve
- 4) ERETIC (digital)

^1H NMR
(60.7 MHz, DMSO- d_6)



qNMR Experiments – Key Considerations




ACCURACY

Weighing
Volumes
Internal calibrant purity



ACQUISITION PARAMETERS

Numbers of points
Acquisition time
Spectral width



RELAXATION

Scan delay must be
chosen to ensure all spins
of interest are fully
relaxed between pulses

Typically, 5-7 times
longest T_1

>> Please see our blog for more details!
Nanalysis website → Resources → Benchtop NMR blog
<https://www.nanalysis.com/nmready-blog>

Development of qNMR Methods

^7Li



Determination of lithium content in brine pools from Chile

^{31}P



Identification of lignin substructures in Canada

^1H



Quantification of MDMA in ecstasy tablets collected at music festivals in the UK

Quantification of Lithium Content in Brine Pools

Problem



- Lithium is becoming increasingly important in our society
- AA & ICP methods require cumbersome sample preps, are difficult to maintain, and are prone to matrix effects

NMR lithium analyzer

Benchtop qNMR for rapid, onsite, simple analysis of lithium content

ROI

- Simple sample preparation
- Reduced operating expenditures
- Rapid & accurate results
- Automated procedure





Quantification of Lithium Content in Brine Pools

Collaboration with Sociedad Chilena de Quimica's (SQM)


16 samples of underground brine pools were obtained

Two sets of these samples were prepared: one was sent to Nanalysis Corp. in Calgary, Canada, while the other remained at SQM in Chile



Method Validation (excerpt)

22-day Long term Stability



44 absolute integrals plotted, (RSD)
(2 data points collected per day for 22 days)

Limit of detection (LOD)
Limit of quantification (LOQ)



LOD (SNR = 3): 40 ppm Li^+
LOQ (SNR = 10): 100 ppm Li^+

Sample Analysis Results


Sample	Li content by NMR	RSD ^a	Li content by AA	Diff (%) ^b
1	0.166	1.86	0.163	1.8
2	0.201	1.06	0.197	2.0
3	0.222	0.83	0.226	-1.8
4	0.292	1.56	0.289	1.0
5	0.349	0.70	0.339	2.9
6	0.467	0.68	0.439	6.4
7	0.478	0.34	0.471	1.5
8	0.617	0.23	0.599	3.0
9	0.684	0.46	0.654	4.6
10	0.811	0.51	0.755	7.4
11	0.850	0.55	0.802	6.0
12	1.136	0.27	1.060	7.2
13	2.110	0.13	1.969	7.2
14	5.271	0.06	5.032	4.7
15	6.380	0.11	6.243	2.2
16	6.480	0.04	6.144	5.5

- Values obtained from NMR are comparable to that of values from AA
- RSD values of replicates were under 2%
- Values differ up to a maximum of 7.2%

Araneda, J. F.; Hui, P.; Riegel, S. D.; Mercado, R.; Green, C. *Analyst* **2021**, 146, 882-888.

Harvey, D. *Analytical Chemistry* 2.1, Open Education Resource (OER) LibreTexts Project Miller, J. N. *Analyst* **1991**, 116, 3-14.

Comparison NMR vs AA



The intercept is 0.0065 (0.0366 and -0.0236)
The slope is 1.0405 (1.029-1.052)
(R^2) is equal to 0.9996

Development of qNMR Methods

^7Li



Determination of lithium content in brine pools from Chile

^{31}P



Identification of lignin substructures in Canada

^1H



Quantification of MDMA in ecstasy tablets collected at music festivals in the UK

Lignin Characterization

Lignin is a by-product of paper production

- Monomers are important building blocks
- Historically obtained from petroleum feedstocks
- Complex cross-linked phenolic polymers
- Structure needs to be determined


Challenge with lignin is structure determination

- Different for each lignin and very complex
- $^1\text{H}/^{13}\text{C}$ NMR analyses suffer from too much overlap


Lignin was successfully derivatized with many nuclei

^{31}P is an excellent nucleus for NMR analysis of lignin

- 100% natural abundance
- Large chemical shift range
- Gives rise to sharp lines



Adapted from Argyropoulos *et al.*



Phosphitylation Approach

This analysis was part of a round-robin study performed to assess a test method described in a future Canadian Standard Assoc. standard method


Phosphitylation of alcohol groups with TMDP

- Cyclohexanol used as internal calibrant

Reaction is immediate and quantitative


- ^{31}P signals in $\delta = 130\text{-}150$ ppm range

Large potential applicability for this approach, but lack of access (for many) to a high-field instrument



^{31}P NMR

(24.3MHz, CDCl_3)



Identification of Lignin Substructures

$^{31}\text{P}\{^1\text{H}\}$ NMR

(Kraft Softwood Lignin, CDCl_3)




Table 1. Chemical shift ranges (^{31}P NMR) used for integration of various regions for the quantification of lignin substructures.

Hardwood and wheat straw	
Functional group	Chemical shift range (ppm)
Aliphatic	150.0–145.4
Internal calibrant	145.4–144.6
β -5	144.5–143.3
4-O-5'	143.3–141.6
5-5'	141.6–140.4
Guaiacyl	140.4–138.5
<i>p</i> -Hydroxyphenyl	138.5–137.0
Carboxylic acid	136.0–133.6
Softwood	
Functional group	Chemical shift range (ppm)
Aliphatic	150.0–145.4
Internal calibrant	145.4–144.6
β -5	144.5–143.0
Syringyl	143.0–141.8
5-5'	141.8–140.4
Guaiacyl	140.4–138.3
<i>p</i> -Hydroxyphenyl	138.3–137.0
Carboxylic acid	136.0–133.6

Identification of Lignin Substructures - Analysis Results

Table 2. Comparison of concentrations (in mmol/g) obtained from 24.3 MHz (Bench) and 162 MHz (HF) ³¹P NMR analysis.

Functional group	KSW			KSW			KHW			OHW			SOD		
	Bench ^a	HF ^b	Error	Bench ^c	HF ^d	Error	Bench	HF	Error	Bench	HF	Error	Bench	HF	Error
X (Aliphatic OH)	1.82 (0.3)	1.71	6%	1.88 (0.7)	1.71	10%	2.02 (1.4)	1.85	9%	2.44 (3.1)	2.37	3%	1.75 (1.1)	1.64	7%
X (5–5')	0.77 (0.4)	0.65	18%	0.72 (1.4)	0.66	9%	0.34 (1.4)	0.33	3%	0.24 (6.1)	0.22	9%	0.34 (3.2)	0.33	3%
X (4-O-5' or Syringyl OH) ^e	0.45 (0.9)	0.41	10%	0.56 (0.5)	0.41	37%	1.90 (2.0)	1.85	3%	1.50 (1.6)	1.51	1%	1.15 (0.5)	1.13	2%
X (β-5)	0.46 (1.8)	0.42	10%	0.50 (1.1)	0.42	19%	0.30 (4.2)	0.26	15%	0.28 (5.6)	0.26	8%	0.28 (1.9)	0.27	4%
X (Guaiacyl OH)	1.80 (1.5)	1.61	12%	1.82 (1.2)	1.63	12%	1.05 (2.2)	1.01	4%	0.80 (1.8)	0.82	2%	1.19 (1.4)	1.15	3%
X (p-Hydroxyphenyl OH)	0.13 (11.9)	0.16	19%	0.15 (3.6)	0.16	6%	0.12 (2.9)	0.19	37%	0.08 (16.1)	0.09	11%	0.46 (5.6)	0.47	2%
X (Carboxylic acid OH)	0.52 (0.7)	0.54	4%	0.56 (2.0)	0.54	4%	0.37 (8.7)	0.41	10%	0.17 (3.7)	0.18	6%	0.98 (0.3)	0.95	3%

Note: The benchtop NMR values are presented as the average of triplicate runs, and the relative standard deviations are provided in parentheses. Results for the individual triplicate analyses are included in the supplementary information (Table S1). ^a1024 scans.

^b128 scans.

^c2048 scans.

^d256 scans.

^e4-O-5' for KHW/OHW and syringyl OH for KSW/SOD. Error (%) = |content by Bench – content by HF|/content by HF x 100.


Development of qNMR methods

^7Li




Determination of lithium content in brine pools from Chile

^{31}P



Identification of lignin substructures in Canada

^1H



Quantification of MDMA in ecstasy tablets collected at music festivals in the UK

Quantification of MDMA in Ecstasy Tablets

Problem

- MDMA quantification important for harm reduction and law enforcement
- Issues of overdoing and deaths related to drug poisoning
- Public Health England: £ 10.7 billion illicit drug costs in 2010/2011 in the UK





Method validation and street sample analysis of 100 tablets collected at music festivals in the UK

- Benchtop qNMR with ethylene carbonate (EC) as internal calibrant

ROI

- Simple sample preparation
- Robust method
- Rapid & accurate results


¹H NMR
(60.1 MHz, DMSO-d₆)



Method Validation (excerpt)

External quality testing

Tablet ID	MDMA content in wt%		Deviation in %
	NMR results	LC-MS results	
16152	47.52	46.0	3.3
18861	31.45	33.3	-5.6
20252	33.95	33.2	2.3
20284	20.06	20.0	0.3
20307LA	30.79	32.7	-5.8
22192/4	38.84	37.1	4.7
24803	2.20	2.25	-2.2
27099	12.92	12.8	0.9
28042	41.53	40.6	2.3
29596	60.69	59.8	1.5




Limit of detection and limit of quantification

- LOD: 0.10 mg/mL (SNR ≥ 3); LOQ: 0.33 mg/mL (SNR ≥ 10)


Street Sample Analysis (excerpt)

Within Batch Variation Testing

Small variation in MDMA content
for the Red Bull tab



High variation in MDMA content
for the Trophy tab



Conclusions

Benchtop technology offers unparalleled NMR access to academic and industrial users by **addressing accessibility issues**

- Low upfront and recurring costs
- No maintenance
- Expert staff not required
- Low infrastructure requirements
- Compact all in one benchtop NMR systems

We have demonstrated highly effective methods using **^7Li** , **^{31}P** , and **^1H qNMR** (for others, please inquire)

- Determination of lithium content in brine pools in Chile
- Identification of substructures in lignin from Canada
- Quantification of MDMA in ecstasy tablets collected at music festivals in the UK



Acknowledgments



Susanne Riegel
Juan Araneda
Matthew Leclerc
Paul Hui
Jonathan Ma
Thais Barbosa
Kiana Nayebi



Simon Rainbow
James Lyons
Thomasine Curzon



Rodrigo Mercado
Christopher Green



Michael Paleologou



Ian Burton



Melanie Pinatton



Anca Frinculescu
Trevor Shine
John Ramsey



Vincenzo Abbate
Nunzianda Frascione

Thank you for your time.

*Please come by our booth
tomorrow and chat with us!*

